

INOGEN ONE® G3

Now that you have unpacked each box, you have everything you need to start using your Inogen One G3 (IO-300). The set-up of your Inogen One G3 is quick and simple. However, before using the Inogen One G3, please read the entire User Manual.



Place your Inogen One G3 in a well ventilated location. Ensure both particle filters are in place.



Insert the Inogen One® G3 battery by sliding battery into place until the latch returns to the upper position.



Connect the AC Power Supply cord to the Inogen One G3.



Plug the Power Supply AC cord into the nearest electrical outlet.



Connect the nasal cannula tubing into the metal nozzle on your Inogen One G3. Oxygen will only be delivered while you are breathing through the nasal cannula connected to the device.



Turn on the Inogen One G3 by pressing the  button.



Set the flow setting to the flow rate prescribed by your clinical professional. The flow setting can be changed by pressing the triangular flow setting buttons that are marked with "+" and "-" symbols. Position the nasal cannula on your face and breathe normally through your nose.



Car/Boat/RV Usage

After you are familiar with operating your Inogen One G3 at home, you are ready to take your Inogen One G3 on the road.

As always, be a careful driver. DO NOT RUSH OR HURRY to set up your Inogen One G3 in the car. Take your time and follow the steps outlined to the right. The Inogen One G3 can operate either by battery or the DC Power Cord while in your car, boat, or RV.



Plug the DC Power cord into the cigarette lighter or auxillary DC power supply.



Plug the power output cord into your Inogen One G3 and position your Inogen One G3 securely inside the car. Operate your Inogen One G3 as normally.

CAUTION

Make sure your Inogen One G3 is secure before driving.



American Medical Sales and Rentals
7032 S Revere Parkway, Suite 320
Centennial, CO 80220
877-774-9271
www.oxygenconcentratorstore.com

Your Equipment Provider:

Your Physician: