

CPAP specifications guide

System One sleep therapy platform

System One	Modes	Comfort features				Ramp (5-45 minutes)	Altitude adjustment
		Flex comfort	Flex demo modes	System One humidity control	System One resistance control		
REMstar	Fixed CPAP	N/A	N/A	Standard	N/A	Adjustable	Manual
REMstar Plus	Fixed CPAP	C-Flex	CPAP C-Flex	Standard	N/A	Adjustable	Manual
REMstar Pro	Fixed CPAP	C-Flex C-Flex+	CPAP C-Flex C-Flex+	Standard	Standard	Adjustable	Automatic
REMstar Pro with AutoIQ	Fixed CPAP Auto-Trial Auto-Check	C-Flex C-Flex+ A-Flex	CPAP C-Flex C-Flex+ A-Flex	Standard	Standard	Adjustable	Automatic
REMstar Auto	Fixed CPAP Auto CPAP	C-Flex C-Flex+ A-Flex	CPAP C-Flex C-Flex+ A-Flex	Standard	Standard	Adjustable	Automatic

System One	Encore patient management						
	Data transfer modes			Data capture			
	Phone-in EncorePro & EncoreAnywhere	SD card EncorePro & EncoreAnywhere	Modem EncoreAnywhere	Patterns of use	Advanced event detection and reporting	Advanced event response	Flow waveforms
REMstar	Hours of use VIC	Standard	Optional	Standard	N/A	N/A	N/A
REMstar Plus	Hours of use VIC	Standard	Optional	Standard	N/A	N/A	N/A
REMstar Pro	Hours of use VIC	Standard	Optional	Standard	AHI, FL, RERA, snore, leak, CA, OA, HY, PB	N/A	Standard
REMstar Pro with AutoIQ	Hours of use VIC	Standard	Optional	Standard	AHI, FL, RERA, snore, leak, CA, OA, HY, PB	FL, RERA, snore, leak, OA, HY	Standard
REMstar Auto	Hours of use VIC	Standard	Optional	Standard	AHI, FL, RERA, snore, leak, CA, OA, HY, PB	FL, RERA, snore, leak, OA, HY	Standard

Legend
 AHI: Apnea Hypopnea Index
 CA: Clear Airway Apnea
 FL: Flow Limitation
 HY: Hypopnea
 OA: Obstructed Airway Apnea

PB: Periodic Breathing
 RERA: Respiratory Effort Related Arousal
 VIC: Visual Inspection Check

PHILIPS

RESPIRONICS

sense and simplicity

Specifications		REMstar DS150S CA151S	REMstar Plus DS250S CA251S	REMstar Pro CA451S	REMstar Pro with AutoIQ DS450S	REMstar Auto DS550S CA551S
Pressure range:	4 to 20 cm H ₂ O	✓	✓	✓	✓	✓
Flex pressure relief:	Mode specific	N/A	✓	✓	✓	✓
Ramp time:	0 to 45 min. (5-min. increments)	✓	✓	✓	✓	✓
Starting ramp pressure:	4 to CPAP; patient-adjustable 4 to min. pressure	✓ N/A	✓ N/A	✓ N/A	✓ (CPAP mode) ✓	✓ (CPAP mode) ✓
Dimensions:	7" L x 5.5" W x 4" H 18 cm L x 14 cm W x 10 cm H	✓	✓	✓	✓	✓
Weight:	3 lbs; 1.36 kg (without humidifier)	✓	✓	✓	✓	✓
Filters:	Foam and optional ultra-fine	✓	✓	✓	✓	✓
Device set-up:	LCD/control wheel/push button	✓	✓	✓	✓	✓
Data storage capacity (minimum):	Display: Compliance VIC 7- and 30-day averages, SD Card: >1 year On-board: up to 6 months	✓	✓	✓	✓	✓
Compliance meter:	Breathing detection	✓	✓	✓	✓	✓
Altitude compensation:	Automatic Manual	N/A ✓	N/A ✓	✓ N/A	✓ N/A	✓ N/A
Electrical requirements:	100 – 240 VAC, 50/60 Hz	✓	✓	✓	✓	✓
Humidification:	Integrates with REMstar System One heated humidifier	✓	✓	✓	✓	✓
Warranty:	Two years; replace with new	✓	✓	✓	✓	✓
Miscellaneous:	Fine pressure calibration Lighted LEDs	✓ ✓	✓ ✓	N/A ✓	N/A ✓	N/A ✓

Ordering information			Ordering information		
Devices	Part number US/CAN	Canada	Accessories	Part number US/CAN	Canada
REMstar with humidifier, SD card	DS150HS	CA151HS	Foam filter kit	1063091	
REMstar with SD card	DS150S	CA151S	Ultra-fine filter kit	1063096	
REMstar Plus with humidifier, SD card	DS250HS	CA251HS	System One heated humidifier	1056209	1056210
REMstar Plus with SD card	DS250S	CA251S	Data management accessories		
REMstar Pro with humidifier, SD card		CA451HS	Part number		
REMstar Pro with SD card		CA451S	SD card (10 pack)	1063859	
REMstar Pro with AutoIQ, with humidifier, SD card	DS450HS		SD card mailer (10 pack)	1065146	
REMstar Pro with AutoIQ, with SD card	DS450S		SD card cover (attached to device)	1063858	
REMstar Auto with humidifier, SD card	DS550HS	CA551HS	SD card reader	1073915	
REMstar Auto with SD card	DS550S	CA551S	Wired flow modem	100200W	
Accessories			Wired flow modem monthly service	100201W	
Part number			Wireless flow modem	100200C	
System One Performance Tube 15 (15 mm)	1074400		Wireless flow modem monthly service	100201C	
Device side panel	1063784				
Water chamber	1063785				
Power supply	1063786				
Carrying case	1063857				
Power cord	1005894				

Respironics, A-Flex, AutoIQ, Auto-Check, Auto-Trial, Encore, EncoreAnywhere, REMstar, and System One are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. The C-Flex mark is used under license.



© 2011 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoech KW 06/29/11 MCI 4104306 PN 1064784

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9666 4444
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics