



# Oxygen simplified

EverFlo concentrators reduce the cost of oxygen therapy.

**PHILIPS**

**RESPIRONICS**

# Ideal oxygen concentrators for you and your patients.

Philips Respironics recognizes the growing need for cost-effective and reliable oxygen concentrators that can reduce the cost of providing oxygen therapy. That's why we've developed a series of forward-thinking oxygen products like EverFlo that simplify your oxygen business.

## Expect industry-leading quality and reliability

EverFlo stationary oxygen concentrators can save you time and money over the course of ownership because of their design and quality. Field-proven over more than two years, EverFlo was rated highest in perceived quality by a 2008 independent survey of homecare providers, and hundreds of thousands of units have been deployed worldwide.\*

Manufactured in the United States, and backed by our deep design knowledge, EverFlo is developed using only high quality parts. With a recessed flow meter, metal cannula connector, strong plastic casing, proven sieve bed, and Pressure Swing Adsorption (PSA) technology, you will notice a reduction in repairs, complaints, and unscheduled visits.

## Enjoy ease of service

EverFlo is designed with fewer parts than other concentrators, resulting in lower overall maintenance costs. EverFlo has just two filters. The external filter can be changed in minutes without special skills or tools and the internal micro disk filter is simply checked during routine maintenance procedures.

## Differentiate your business

EverFlo concentrators are backed by value-added tools. From in-person and online training to financing options, you have what you need to succeed. And, since Philips Respironics offers a wide range of sleep and home respiratory products and comprehensive support, we can be the source for many of your product needs.

\* HME Home Oxygen Survey, 2008. Data on file. Philips Respironics.



## Easy to store, transport, and deliver.

EverFlo is compact, has a streamlined and integrated handle, and weighs only 31 pounds. This results in lower shipping and storage costs and reduced risk of worker injury.



## Give your patients what they want.

Your patients will appreciate the compact non-medical design of the EverFlo oxygen concentrator, because it will fit right into their home environment. The unit is designed to be small, light, easy to transport and will fit into the trunk of a car. EverFlo has an easy-to-mount humidifier bottle connection and no patient filter.



Molded handle that's easy to pick up and transport with one or both hands



Humidifier bottle holder compatible with all bottle styles

## EverFlo product specifications

|  |  |
|--|--|
| Input voltage                                  | 120 VAC +/- 10%  |
| Input frequency                                | 60hz   |
| Average power consumption                      | 350 watts  |
| Oxygen concentration* (at 5 lpm)               | 93% (+/- 3%)   |
| Liter flow                                     | 0.5-5 liters per minute                                    |
| Weight   | 31 lbs (14 kg)   |
| Sound level                                    | 45 dBA typical   |
| Dimensions                                     | 23" (584 mm) h, 15" (381 mm) w, 9.5" (241 mm) d            |
| Required maintenance                           | Compressor intake filter: change every 2 years             |
| OPI (oxygen percentage indicator) alarm levels | Low oxygen: 82%, very low oxygen: 70%                      |
| Operating temperature                          | 55°F to 90°F (12°C to 32°C)                                |
| Storage/transport humidity                     | -30°F to 160°F (-34°C to 71°C) up to 95% relative humidity |
| Operating humidity                             | Up to 95% relative humidity                                |
| Outlet pressure                                | 5.5 PSI  |
| Operating altitude                             | 0 To 7,500' (0 to 2,286 m)                                 |

## EverFlo ordering information

| Concentrator systems   | Part number |
|--|-------------|
| EverFlo oxygen concentrator  | 1020000     |
| EverFlo oxygen concentrator with OPI (oxygen percentage indicator)   | 1020001     |
| EverFlo Q oxygen concentrator  | 1020015     |
| EverFlo Q oxygen concentrator with OPI (oxygen percentage indicator) | 1020014     |
| Concentrator accessories   | Part number |
| Compressor intake filter   | 1038831     |
| Micro disk filter  | 520         |
| Humidifier connector tube  | 1039642     |
| Stand-alone pediatric flow meter                                     | H537-1      |
| Integrated pediatric flow meter                                      | H644        |
| Flow meter locking kit   | 528         |
| EverFlo manual in English, French, Spanish, and Portuguese           | 1020121     |
| EverFlo service manual   | 1038571     |
| Extended warranty  | Part number |
| Two-year extended warranty with compressor                           | 1038571     |

\* Device operation above or outside of the specified voltage, LPM, temperature, humidity, and/or altitude values may decrease oxygen concentration levels.



Internal tubing storage keeps tubing conveniently available for use



Metal cannula connector for greater durability



Interface panel with on/off switch and alarm indicators



Recessed flow meter reduces the risk of breakage

# Freedom of choice

We offer a full line of oxygen concentrators.

## **EverFlo Q ultra quiet oxygen concentrator with all of the features of our EverFlo, plus:**

- Provides ultra-quiet premium sound level of 40 dBA
- Reduces noise complaints and associated costs and time
- Differentiates your oxygen business to physicians

## **EverGo portable oxygen concentrator (shown below):**

- Eight-hour battery life
- 1050 ml/min oxygen capacity
- Less than ten pounds
- Easy-to-use non-medical design that plugs into any AC or DC power source
- FAA approved for air travel and permitted for in-flight use



Pair EverFlo with EverGo for an ideal oxygen-generating solution.

EverFlo provides a low-maintenance, low-cost concentrator for stationary use and EverGo is a lightweight portable concentrator with eight-hour battery life. The EverFlo and EverGo system can reduce your cost of providing oxygen therapy by:

- Eliminating recurring delivery costs and cylinder hassles
- Reducing maintenance costs
- Easing inventory management
- Minimizing patient calls

**Philips Healthcare is part of  
Royal Philips Electronics**

**How to reach us**

[www.philips.com/healthcare](http://www.philips.com/healthcare)

[healthcare@philips.com](mailto:healthcare@philips.com)

fax: +31 40 27 64 887

**Asia**

+852 2821 5888

**Europe, Middle East, Africa**

+49 7031 463 2254

**Latin America**

+55 11 2125 0744

**North America**

+1 425 487 7000

800 285 5585 (toll free, US only)

Philips Respironics  
1010 Murry Ridge Lane  
Murrysville, PA 15668

Customer Service  
+1 724 387 4000  
800 345 6443 (toll free, US only)

Respironics Asia Pacific  
+65 6298 1088

Respironics Australia  
+61 (2) 9666 4444

Respironics Deutschland  
+49 8152 93 06 0

Respironics Europe, Middle East, Africa  
+33 1 47 52 30 00

Respironics France  
+33 2 51 89 36 00

Respironics Italy  
+39 03 62 63 43 1

Respironics Sweden  
+46 8 120 45 900

Respironics Switzerland  
+41 6 27 45 17 50

Respironics United Kingdom  
+44 800 1300 845

[www.philips.com/respironics](http://www.philips.com/respironics)

Respironics, EverFlo, and EverGo are trademarks of Respironics, Inc. and its affiliates.



©2009 Koninklijke Philips Electronics N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.  
CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Geyer WMB 06/08/09 MCI 4102133 PN 1038667