


# FREESTYLE® COMFORT®

Portable Oxygen Concentrator



## The future is here

With its ergonomic shape, the lightweight and affordable FreeStyle® Comfort® rests against the natural curves of your body, making it comfortable to wear on the go. FAA-approved for commercial air flights and equipped with user-replaceable batteries, this new oxygen therapy device allows you to extend your time away from home. Our optional telehealth solution — CAIREview™ powered by SynsorMed\* — allows for monitoring of the device's location and performance, as well as remote troubleshooting, all from the convenience of a desktop computer or smart device.



### A LIGHTWEIGHT SOLUTION

Weighing only 5 lbs., the FreeStyle® Comfort® can go anywhere you want to go.


Enjoy the convenience of wearing it with or without your carry bag by simply attaching the no-slip, padded shoulder strap.



### BREATHE EASY, WORRY LESS

UltraSense™, our sensitive breath detection technology, ensures oxygen is delivered to you when you need it.

Bluetooth technology connects you to CAIREview™ — a telehealth solution that your provider can use to help ensure you get the support you need.



### COMFORTABLE TO WEAR, EASY TO USE

The ergonomic, curved silhouette rests against the natural curves of your body, making it more comfortable to wear when you are out and about.

User-replaceable batteries allow you to extend your time away from home and stay active with up to 4 hours of battery life.

\*CAIREview™ is not intended for active, real-time patient monitoring.

# FREESTYLE® COMFORT®

Portable Oxygen Concentrator

## Includes

- Portable Concentrator
- Carry Bag
- Battery Pack
- AC/DC Universal Power Supply and Cords
- Shoulder Strap
- UltraSense™ technology
- Other SKUs may include extra accessories

## Product Features

- User-replaceable battery pack with up to 4 hours of battery life at setting 2
- Easy to wear with a shoulder strap
- FAA-approved for commercial airline use
- Fully functional on wall (AC) or motor vehicle (DC) power
- Easy-to-use full-color interface
- Can come equipped with the optional telehealth solution



## Specifications

**FLOW RATE**  
Pulse 1-5

**DIMENSIONS**  
10.0 x 7.3 x 3.1 in  
(25.4 x 18.5 x 7.9 cm)  
8-cell battery pack


11.0 x 7.3 x 3.1 in  
(27.9 x 18.5 x 7.9 cm)  
16-cell battery pack

**WEIGHT**  
5 lbs (2.3 kg)  
8-cell battery pack  
6 lbs (2.7 kg)  
16-cell battery pack

**SOUND LEVEL**  
43 dB(A) @ Setting 2

**BATTERY DURATION**  
8-cell battery at setting  
2 – up to 4 hours  
16-cell battery at setting  
2 – up to 8 hours

**WARRANTY**  
3 years FreeStyle® Comfort®  
(Includes sieve beds)  
1 year battery and power supplies



## Available Accessories



**8-CELL BATTERY PACK**



**16-CELL BATTERY PACK**



**CARRYING BAG**



**DESKTOP BATTERY PACK CHARGER (OPTIONAL)**



**AC POWER SUPPLY (ALL CORDS INCLUDED)**



**DC POWER CORD**



**EXTERNAL GROSS PARTICLE FILTERS**

### OTHER ITEMS AVAILABLE:

- Shoulder Strap
- Carry-all Accessory Bag

### American Medical Sales and Rentals

7032 S Revere Parkway, Suite 320  
Centennial, CO 80112  
Phone: 877-774-9271

[www.oxygenconcentratorstore.com](http://www.oxygenconcentratorstore.com)  
E-mail: [sales@amsrco.com](mailto:sales@amsrco.com)

