

# AirSep® FreeStyle™

PORTABLE OXYGEN CONCENTRATOR

## FREEDOM HAS A NEW LOOK

CAIRE introduces updates to the original lightweight portable oxygen concentrator — the AirSep® FreeStyle™. The new design adds features and benefits that will not only make the device more appealing to patients, but even more reliable for providers.

- **New Power Supply Connection**

Completely re-designed connector featuring **multi-directional installation** with **no pins to align** and **no release button to press** is simpler than ever to use. **Robust metal connection points** ensure the highest durability of both the power supply and FreeStyle for patients and provider alike.

- **New User-Replaceable Battery**

**Swap battery packs on the go** instead of waiting to re-charge an internal battery or having to add heavy external batteries. The battery pack will **re-charge anytime AC or DC power is connected** and is easily accessible from the front of the device.

- **Longer Battery Life**

The new battery pack provides an additional **1 hour of duration** at setting 2, providing up to **3 ½ hours with a single battery pack at setting 2.**

- **New Color Scheme**

New color scheme is **more appealing** and **non-medical** in appearance. The matte finish will help **resist everyday scratches.**



  
**CAIRE**®  
A Chart Industries Company

 **OXYGENConcentratorStore**  
FROM AMERICAN MEDICAL SALES & RENTALS

Toll Free 1.877.774.9271

[www.OxygenConcentratorStore.com](http://www.OxygenConcentratorStore.com)

  
**CHART**®

## PROVIDER FEATURES

- Settings 1-3 cover majority of patients without the added weight or cost of larger capacity units.
- Non-delivery modality that eliminates expensive deliveries.
- Time-proven and trusted quality
- Standard 3-year warranty includes sieve beds
- Increase referrals and prescriptions

## PATIENT FEATURES

- Under 5 pounds with a small compact frame.
- User-replaceable battery pack with up to 3½ hour battery life at setting 2.
- Easily wearable with a handle, shoulder strap, or backpack.
- FAA Approved
- Fully functional on AC or DC power.
- Simple, One-Touch Operation

## CLINICAL FEATURES

- UltraSense™ Technology is the most sensitive breath detection on the market.
- Quickly detects and pulses oxygen to ensure delivery in the first 1/3 of the breath.
- Adjusts bolus size based on breath rate.
- Lightweight package promotes greater ambulation and quality of life.



## SPECIFICATIONS

<b>Flow Rates</b>	Pulse 1-3
<b>Dimensions</b>	8.6 x 6.1 x 3.6 in (21.8 x 15.5 x 9.1 cm)
<b>Weight (w/Battery)</b>	4.9 lbs (2.2 kg)
<b>AC Power</b>	100-240 VAC, 50-60Hz (1.75 amp)
<b>DC Power</b>	11-16 VDC (8.0 amp max)
<b>Operating Temp.</b>	41°F to 104°F (5°C to 40°C)
<b>Sound Level</b>	41 dB(A) @ Setting 2
<b>Operating Altitude</b>	Up to 12,000 ft (3657.6 m)
<b>O2 Concentration</b>	90% (+5.5% / -3%)
<b>Battery Duration</b>	Setting 1 - 6.5 hours (6:30) Setting 2 - 3.5 hours (3:30) Setting 3 - 3.0 hours (3:00)
<b>Battery Re-Charge</b>	4 - 5 hours
<b>Warm-Up Time</b>	< 2 minutes
<b>O2 Output</b>	Setting 1 - 227.5 mL/min Setting 2 - 420 mL/min Setting 3 - 500 mL/min
<b>Limited Warranty</b>	3 years unit 1 year battery & power supplies

## ORDERING INFORMATION

**FreeStyle Includes:** Portable Concentrator, Battery Pack, AC/DC Universal Power Supply, AC Power Cord, DC Power Cable, Carrying Case, Shoulder Strap, User Manual, Carry-all Accessory Bag

<b>AS095-101</b>	USA
<b>AS095-102</b>	European
<b>AS095-B101</b>	USA - 2 batteries
<b>AS095-B102</b>	Euro - 2 batteries
<b>AS095-103</b>	Australia
<b>AS095-104</b>	UK
<b>AS095-B104</b>	UK - 2 batteries
<b>AS095-105</b>	Brazil
<b>AS095-106</b>	Swiss

*\*other countries' power supplies may be available - contact customer service\**

## ACCESSORIES

<b>BT033-1</b>	Spare Battery Pack
<b>MI284-1</b>	Backpack Harness
<b>MI419-1</b>	Carrying Case
<b>MI320-1</b>	Carry-All Bag
<b>MI426-1</b>	Shoulder Strap
<b>F1144-1</b>	Intake Filter
<b>PW032-1</b>	Universal Power Supply (US)
<b>PW032-2</b>	Universal Power Supply (EU)
<b>PW032-4</b>	Universal Power Supply (UK)
<i>Power Supply Included</i>	
<b>BT032-1</b>	Desktop Charger (US)
<b>BT032-2</b>	Desktop Charger (EU)
<b>BT032-4</b>	Desktop Charger (UK)
<i>Cords Included</i>	

  
**CAIRE**®  
A Chart Industries Company

 **OXYGENConcentratorStore**  
FROM AMERICAN MEDICAL SALES & RENTALS

Toll Free 1.877.774.9271

[www.OxygenConcentratorStore.com](http://www.OxygenConcentratorStore.com)

  
**CHART**®