



OXYGENCONCENTRATORSTORE

FROM AMERICAN MEDICAL SALES & RENTALS

Home Oxygen Concentrator Comparison

Unit Name	Oxygen Flow Options	Continuous Flow Liters	Weight	Dimensions (WxHxD) inches
Caire Companion 5 	Continuous	0.5 to 5 LPM	36 lbs (or 16.3 kg)	12.5 x 21.5 x 13.5 in
Respironics EverFlo 	Continuous	0.5 to 5 LPM	31 lbs (or 14.1 kg)	15 x 22.8 x 9.5 in
OxyGo OxyHome 	Continuous	0.5 to 5 LPM	41.8 lbs (or 19 kg)	13.4 x 25.6 x 11.8
Inogen at Home 5L 	Continuous	1 to 5 LPM	18 lbs (or 8.2 kg)	13 x 16.5 x 7 in
Respironics Millennium M10 	Continuous	1 to 10 LPM	53 lbs (or 24 kg)	13.3 x 26.8 x 18.9 in
Caire Intensity 10 	Continuous	2 to 10 LPM	58 lbs (or 26.3 kg)	16.5 x 27.5 x 14.5 in



Unit Name	Noise	Humidifier Bottle Compatible	CPAP Machine Compatible	Power Consumption	Max Operating Altitude	Warranty
Caire Companion 5 	45 dBA at 2 LPM	Yes	Yes	Less than 350 W	1,253 to 9879 feet	Three-Year Manufacturer's Warranty
Respironics EverFlo 	43 dBA	Yes	Yes	Less than 350 W	0 to 7,500 feet	Three-Year Manufacturer's Warranty
OxyGo OxyHome 	Less than 45 dBA	Yes	Yes	Less than 350 W	0 to 10,000 feet	Three-Year Manufacturer's Warranty
Inogen at Home 5L 	Less than 50 dBA	Yes	Yes	Less than 275 W	0 to 10,000 feet	Three-Year Manufacturer's Warranty
Respironics Millennium M10 	Less than 50 dBA	Yes	Yes	Less than 600 W	0 to 1,368 feet	Three-Year Manufacturer's Warranty
Caire Intensity 10 	58 dBA	Yes	Yes	Less than 600 W	-1,250 to 10,000 feet	Three-Year Manufacturer's Warranty