

So much more than just  
a medical device.

**OxyHome®** allows you  
to enjoy life at home—  
without interruption.



**American Medical Sales and Rentals**

15055 E. Hinsdale Drive, Bldg D  
Centennial, CO 80112  
Phone: +1 (877) 774-9271

[www.oxygenconcentratorstore.com](http://www.oxygenconcentratorstore.com)  
E-mail: [sales@amsrco.com](mailto:sales@amsrco.com)



Connect With Us



[www.oxygenconcentratorstore.com](http://www.oxygenconcentratorstore.com)

**Oxy  
HOME®**

5L Stationary Concentrator

So quiet... you won't  
even know it's on.



PN B150-1023

[www.oxygenconcentratorstore.com](http://www.oxygenconcentratorstore.com)

The Quiet Concentrator



## The Future of Stationary Oxygen

The new OxyHome® Stationary Concentrator seamlessly blends style and functionality with its innovative design. Its compact size ensures that it fits effortlessly into any living space. But what truly sets it apart is its near-silent operation.

Enjoy daily activities without interruption; OxyHome is so quiet, you won't even notice it's running!

### OxyHome Benefits

- Quiet operation
- Small footprint
- Purity check by percentage
- Alarm indicators:  
no flow, low/high pressure, purity, power
- Twin blower cooling design
- Contemporary look
- 3-year replacement warranty



### OxyHome Specifications

Dimensions	13.4"(W) x 11.8"(D) x 25.6"(H) 34cm x 30cm x 65cm
Weight	41.8 lbs (19 kg)
Flow Spec.	0 - 5 LPM
Oxygen Concentration	93% (+3%, -3%)
Oxygen Output Pressure	8 (+2, -2) psi
Rated Voltage/Frequency	120 VAC ± 10% 60Hz
Power Consumption	350 Watts
Sound	≤ 45 dB(A)

- ✓ Enjoy uninterrupted, peaceful sleep without disturbances from a noisy machine
- ✓ Move freely around your home, engaging in daily activities
- ✓ Interact with family and friends without the constant background noise from your oxygen concentrator