

# Oxy HOME<sup>®</sup>

## Stationary Concentrator

Keep Going.<sup>sm</sup>

NEW



### NEW BENEFITS + IMPROVEMENTS!

- **Auxiliary Oxygen Port** – Compatible with industry standard home fill units
- **Simplified User Control Panel**
- **Compatible with External Pediatric Flow Meters**
- **New Push-Style Power Button**
- Quiet Operation and Small Footprint
- Alarm Indicators: No Flow, Low Pressure, Purity, Power
- Twin Blower Cooling Design
- Contemporary Design
- 3-Year Replacement Warranty

### PRODUCT SPECS Part Number 1400-5100

|                                |  |
|--------------------------------|--|
| <b>Flow Specifications</b>     | 0 - 5 LPM  |
| <b>Oxygen Concentration</b>    | 93% (+3%, -3%)                                       |
| <b>Oxygen Output Pressure</b>  | 8 (+2, -2) psi                                       |
| <b>Rated Voltage/Frequency</b> | 120 VAC ±10% 60Hz                                    |
| <b>Power Consumption</b>       | 350 Watts  |
| <b>Sound</b>                   | ≤ 45 dB(A)   |
| <b>Dimensions</b>              | 13.4" W x 11.8" D x 25.6" H<br>34 cm x 30 cm x 65 cm |
| <b>Weight</b>                  | 41.8 lbs (19 kg)                                     |



CONTACT US TODAY TO LEARN MORE

1 (877) 774-9271 [www.oxygenconcentratorstore.com](http://www.oxygenconcentratorstore.com)

PN M538-0125