


# Inogen Rove 6



Inogen Rove 6  
AC Power Supply\*  
Item #RP-109



Inogen Rove 6  
DC Power Cable\*  
Item #BA-306



Inogen Rove 6  
Standard Battery\*  
Up to 6:15 hours run time  
Item #BA-500



Inogen Rove 6  
Carry Bag\*  
Item #CA-500



**FLOW  
SETTINGS  
1-6!**

## The Inogen Rove 6 System (IS-501):

- Compact and Lightweight
- Flow Settings 1-6
- Intelligent Pulse Oxygen Delivery
- Designed for 24 Hour Use
- Easy-to-read LCD display

\*Included with system

\*\*Accessory sold separately



Inogen Rove 6  
Extended Battery\*\*  
Up to 12:45 hours run time  
Item #BA-516



Inogen Rove 6  
Backpack\*\*  
Item #CA-550

# Inogen Rove 6



## Introducing the Inogen Rove 6: Your Oxygen Companion

Experience oxygen liberation with the Inogen Rove 6, the pinnacle of portable oxygen delivery innovation. Weighing in as the most oxygen-efficient portable concentrator available, the Inogen Rove 6 offers unparalleled performance in a compact design. With customizable pulse flow settings ranging from 1 to 6, it empowers users with newfound independence around the clock.

Whether you're at home or on the move, the Inogen Rove 6 keeps pace with your lifestyle, offering the freedom of Oxygen. Say goodbye to multiple solutions – with the Inogen Rove 6, you have everything you need for oxygen therapy, whether you're at home or away.

<b>Weight</b>	4.8 lbs. (includes standard battery)
<b>Size</b>	Length: 7.18" Width: 3.27" Height: 8.14" (includes standard battery)
<b>Oxygen Flow</b>	Pulse dose delivery system Six flow settings 1, 2, 3,4,5,6
<b>Power</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: 13.5-15.0 VDC, 100W (for mobile use in car)
<b>Battery</b>	Duration (standard battery): Up to 6 hours 15 minutes Recharge up to 3 hours with AC or DC power Duration (extended battery): Up to 12 hours 45 minutes Recharge up to 4 hours with AC or DC power
<b>Noise Level</b>	39 dBA* typical at setting 2
<b>Warranty</b>	3 year manufacturer warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display Inogen Connect App
<b>Use</b>	Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft

\*At flow setting 2

\*\*See Inogen Warranty Statement for complete warranty information

### American Medical Sales and Rentals

15055 E. Hinsdale Drive, Bldg D  
Centennial, CO 80112  
Phone: +1 (877) 774-9271

[www.oxygenconcentratorstore.com](http://www.oxygenconcentratorstore.com)  
E-mail: [sales@amsrco.com](mailto:sales@amsrco.com)

