



ResMed

Mobi accessories



Mobi external battery

Offering up to 4.5 hours of additional running time at setting 1, the Mobi external battery allows you to get out and enjoy life longer.



Mobi backpack

The Mobi backpack keeps your machine, batteries and nasal cannula secure and out of sight, while the backpack straps conveniently conceal the cannula.



Mobi mounting straps (4)

A set of four adjustable straps feature hook and clip attachments to help you securely attach Mobi to a scooter, walker or wheelchair.



Nasal cannula

7-ft silicone nasal cannula offers a comfortable fit and ample length for freedom of movement.

*Accessories sold separately.

Technical specifications

Mobi device

Weight	5.5 lbs (2.49 kg)		
Size (W x H x D)	9.60" x 7.95" x 4.49" (244 mm x 202 mm x 114 mm)		
Oxygen delivery	Pulse oxygen therapy		
Device settings	1, 2, 3, 4		
Battery life*	Setting	Internal battery	External
	1	8.6 hrs	4.5 hrs
	2	6 hrs	3 hrs
	3	3.8 hrs	1.9 hrs
	4	3.2 hrs	1.6 hrs

*Hours are approximate and based on a breath rate of 20 bpm. The device will operate up to the battery life (hours) specified.

ResMed.com/Mobi

Distributed by ResMed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA +1 858 836 5000 or 1 800 424 0737 (toll free). See ResMed.com for other ResMed locations worldwide. Mobi and Pulse-Wave are trademarks and/or registered trademarks of the ResMed family of companies. Specifications may change without notice. For patent and other intellectual property information, see ResMed.com/ip. © 2018 ResMed. 10110452/Z 2018-05



ResMed

Mobi™



1-877-774-9271 | www.oxygenconcentratorstore.com

Get the most out of life

ResMed Mobi is a portable oxygen concentrator for everyday living. It features a lightweight design, long battery life and delivers the oxygen you need to live your life.



Go further

Mobi is designed to help you enjoy an active lifestyle. From daily errands to long-distance travel, Mobi can go where you go.

Oxygen for you

Mobi uses Pulse-Wave™ technology to efficiently deliver your selected dose of oxygen each time you inhale. As your activity level changes, Mobi senses your breathing and automatically switches between active and rest mode.

Simple to use

Mobi is designed for one-touch startup and simple setting changes. An easy-to-read display shows battery life and indicates when each dose of oxygen is delivered.

Casual style

The Mobi carry bag is both stylish and functional. It features a casual but sophisticated design that conceals tubing and batteries. Or choose the Mobi backpack – great for day trips and travel.

A name you can trust

ResMed has been changing lives with every breath for nearly thirty years. Mobi is the latest addition to our award-winning family of respiratory solutions.



Mobi™
Portable oxygen concentrator

*Accessories sold separately.

1-877-774-9271 | www.oxygenconcentratorstore.com