

> **HOW DO I GET ONE?**

- 1. Tell your doctor about the LifeChoice Activox POC model you are interested in. If the POC will meet your oxygen needs, ask that the prescription be written specifically for the "LifeChoice Activox Portable Oxygen Concentrator" with the direction "Do Not Substitute."
- 2. Contact your local provider, or find one in your area by calling 1.800.220.0977 or by visiting [www.InovaLabs.com](http://www.InovaLabs.com).



**PORTABLE OXYGEN THAT'S  
READY TO GO WHEN YOU ARE**

**1.877.774.9271**  
[www.oxygenconcentratorstore.com](http://www.oxygenconcentratorstore.com)



**NO TANKS  
NO HASSLES  
NO BOUNDARIES**

> **STAY ACTIVE**

LifeChoice Activox POCs are lightweight, reliable and ready to go—for life with no compromises. Designed to be ready when you are—there is no waiting for tanks to fill and no heavy carts to drag around.

> **BEST-IN-CLASS PORTABILITY**

Available in three lightweight models with the longest-lasting internal batteries available, LifeChoice Activox POCs are the ideal portable oxygen therapy for active lifestyles.



Feature	Model		
	4L	Pro	Sport
Weight	4.8 lbs.	4.3 lbs.	3.9 lbs.
Oxygen Delivery	PULSE-WAVE™ Delivery		
	1-4 LPMeq	1-3 LPMeq	1-3 LPMeq
Internal Battery Runtime	Up to 10.25 Hours	Up to 12.5 Hours	Up to 5.5 Hours
FAA Approved	✓	✓	✓
Tank-Free	✓	✓	✓



> **LIGHTWEIGHT AND HANDS-FREE**

LifeChoice Activox POCs are lightweight and do not require the use of a cart. Plus, each model comes equipped with a convertible carry case that allows you to wear the system as a backpack, shoulder bag, waist pack or briefcase. These options empower you to be hands-free so you can complete activities with comfort and ease.

> **PERFECT FOR TRAVEL**

All models (4L, Pro and Sport) are FAA-approved for domestic and international airline travel. With unmatched battery life and unique delivery technology to keep you saturated whether your home or away, LifeChoice Activox POCs are the perfect travel companion for active oxygen users.



> **PEACE OF MIND OPERATION**

As a tank-free solution, you never have to worry about running out of oxygen. Recharge the battery at home or on the go—even while you’re using it! Simple controls adjust the settings while clearly stated alarm and charging notifications keep you informed during operation. Each model also features an advanced oxygen sensor to let you know when purity levels change.

The 4L Model features an advanced transfective display, as shown below, for easy viewing in bright sunlight.



> **INNOVATIVE AUTO MODE & PULSE-WAVE™ TECHNOLOGY**



Every LifeChoice Activox POC features powerful Auto Mode and PULSE-WAVE™ technologies which work together to keep you saturated around the clock.

Auto Mode detects changes to your breath rate and automatically adjusts the delivery of oxygen, while PULSE-WAVE™ Technology delivers that oxygen to you in a manner that is congruent with your breath rate so that all of the oxygen dose is received. This type of delivery is gentle and helps minimize uncomfortable side effects associated with other forms of oxygen delivery such as dry nose and nose bleeds.



> **EVERYTHING YOU NEED IN ONE AFFORDABLE PACKAGE**

The LifeChoice Activox POC is reimbursed by Medicare and most private insurance plans. The Standard Package includes all the necessary elements for easy use and portability.



- 1 LifeChoice Activox 4L Model Shown**

**2 4-Way Carry Case**  
Use as a backpack, shoulder bag, waist pack or briefcase

**3 Adjustable Straps**  
For use with 4-Way Carry Case

**4 DC Power Supply**  
For charging on the go in vehicle
- 5 AC Power Supply**  
For charging at home or on the go through a standard wall outlet

**6 Standard 7-Foot Long Single Lumen Nasal Cannula\***

**7 Accessory Bag**  
To hold power supplies, extra straps and manual



**Additional External Battery Available†**

**Standard 3-Year Limited Warranty**  
(One year warranty for Sieve Bed, batteries and accessories)

\*Nasal cannula may not be included.  
†Not included in Standard Package. Sold separately.

> **WHICH MODEL IS RIGHT FOR YOU?**

Feature	Model		
	4L	Pro	Sport
Weight	4.8 lbs.	4.3 lbs.	3.9 lbs.
Oxygen Purity	90% ± 3%		
Oxygen Delivery	PULSE-WAVE™ Delivery		
	I-4 LPMeq	I-3 LPMeq	I-3 LPMeq
Internal Battery Runtime*	10.25 H @ 1 8.25 H @ 2 5 H @ 3 4 H @ 4	12.5 H @ 1 7.25 H @ 2 6.25 H @ 3	5.5 H @ 1 4 H @ 2 3.25 H @ 3
External Battery Runtime*†	4.75 H @ 1 3.75 H @ 2 2.75 H @ 3 2.25 H @ 4	3 H @ 1 3 H @ 2 2 H @ 3	3 H @ 1 3 H @ 2 2 H @ 3
Charge Time*	4 H Internal 2 H External	4.5 H Internal 3 H External	3 H Internal 3 H External
Average Sound @ 2	44 dB	46 dB	46 dB
Size	9.05" w x 7.875" h x 4.38" d		
Display Type	Transfective w/Buttons	Standard Buttons	Standard Buttons
Warranty	Standard 3-Year Limited		
FAA Approved	✓	✓	✓

\*Hours are approximate and based on breaths per minute.  
†Not included in Standard Package.