



Accessories that support an active lifestyle for ambulation & mobility, travel & more.

We offer our customers a turnkey solution for ambulation!

power cartridge (battery)

Item# 7082



- New, easy-to-use, black button features intuitive downward release
- Provides up to 5.0 hours of power
- Charges using AC/DC power, or the Eclipse Desktop Charger (#7112)
- Comes standard with each Eclipse
- **File HCPCS code E1356**

AMBULATION & MOBILITY

desktop charger

Item #7112



- Charges the Eclipse Battery (#7082) in less than 3 hours
- Use anywhere in the world, by combining with the Universal Power Plug Set (#5591)
- Universal Power Plug Set sold separately
- **File HCPCS code E1357**

eclipse universal cart

Item #5991



- Features six-inch wheels for higher ground clearance
- Ideal for all seasons, all terrain
- Easy access to Battery
- Handle slides up and down with different settings for personalized adjustment
- **File HCPCS code E1354**

smartpulse oximeter

Item# 5468



- Four-year warranty
- Multi-colored LED display
- Audio alert when in "POOR" range
- FDA approved & made in the USA
- Includes neck strap, carrying case and two AAA batteries

eclipse 3 AC adapter with NEMA power cord

Item #5941



- Intended for use at home, or wherever standard AC power is available
- SeQual branded to help distinguish from E1 and E2 power supplies
- Green LED illuminates when AC Power Supply is supplying 28VDC power

eclipse wheelchair pack

Item #5220



- Attaches to any push-handle wheelchair using adjustable, buckle-style clips
- Eclipse fits easily into the main compartment
- Velcro flap provides easy access to the Eclipse Battery (#7082)



TRAVEL

travel accessory kit

Item #5093



- SeQual's featured bundle accessory
- Contents include Eclipse Travel Case (#5010), Eclipse Battery (#7082), Desktop Charger (#7112) and North American Power Cord (#3890)
- File HCPCS codes E1356 (battery) and E1357 (charger)

travel case

Item #5010



- Features adjustable Velcro® inner partitions
- Neatly stores Eclipse accessories
- Rubber handle offers a sturdy grip
- Easily slips over the handle of all Eclipse Carts

eclipse protective cover

Item #5052



- Reduces external wear and tear
- Protects against spillage and inclement weather
- Secures with convenient zipper closure
- Clear plastic window allows access to control panel

universal power plug set

Item #5591



- Ideal for international travel with the Eclipse
- Europe, UK, New Zealand, North/South America, Japan and Taiwan
- Rated at 10A to support the current demand of the Eclipse

eclipse 3 DC adapter

Item #5942



- Intended for use with DC accessory outlets, such as those found in motor vehicles
- Green LED illuminates when DC Power Supply is supplying 26VDC power
- File HCPCS code E1358

eclipse accessory bag

Item #7104



- Neatly stores Eclipse accessories
- Rubber handle offers a sturdy grip
- Easily slips over the handle of all Eclipse Carts

AUXILIARY & MORE

eclipse data acquisition tool

Item# 5535



- Enables firmware upgrades on-site
- Reduces downtime and maintenance costs
- Triage the device for equipment fault vs. user error
- Print and retain historical event logs

Pelican® hard shell courier case

Item #4204



- Protects contents during repeated transport
- Watertight and shock absorbent
- Includes the Eclipse Travel Case (#5010)
- Lifetime warranty

humidifier adapter kit

Item #7116



- Attaches humidifier bottle to Eclipse unit (humidifier bottle not included)

1-877-774-9271

www.oxygenconcentratorstore.com

