

American Medical Sales and Repair has been serving customers requiring oxygen needs for more than a decade. We offer the lowest prices and the best customer service. Our office is open 7 days a week to assist you!

Renting a portable oxygen concentrator can provide an oxygen patient with freedom and mobility.

- All units run on AC power, DC power (automobile cigarette lighters), and batteries
- All units are FAA approved
- Easy rental process

**PRICING:**

- Portable Concentrators with AC/DC chargers & 1 battery start at \$245 per week (plus shipping)
- Additional batteries and accessories are an extra charge

**DELIVERY:**

- Immediate delivery available in the Metro Denver area — 7 am - 9 pm
- UPS Ground, 2nd Day Air and Overnight Shipping. (Expedited Delivery not available in all areas)

**AMSR ALSO OFFERS:**

- New and used units for sale
- An attractive Rent-to-Own program
- Home Units starting at \$15 a day

**TYPES OF UNITS WE RENT:**

			
Respironics EverGo	Respironics SimplyGo	Inogen One G2	SeQual Eclipse 3
10 pounds	10 pounds	8 pounds	18 pounds
Pulse Delivery Only (settings 1 - 6)	Continuous Flow (.5 - 2 LPM)	Pulse Delivery Only (settings 1 - 5)	Continuous Flow (.5 - 3 LPM)
	Pulse Dose Delivery (settings 1 - 6)		Pulse Dose Delivery (settings 1 - 6)

We also rent the Inogen One G3, SeQual Eclipse 5, InovaLabs Lifechoice and Oxlife Independence.

**BATTERY LIFE (IN HOURS):**

Battery life varies with the setting. Use this convenient table to estimate the number of batteries you will need.

Pulse Dose Machines	Setting 1	Setting 2	Setting 3	Setting 4	Setting 5	Setting 6
Respironics Evergo (1 battery)	6	4	3	2	2	2
Inogen One G2 - 12 cell	4	3.5	3	2.5	2	NA
Inogen One G2 - 24 cell	8	7	6	5	4	NA
SeQual Eclipse 3 (pulse)	4.4	3.6	3	2.6	2.3	2.1
Respironics SimplyGo (pulse)	3.4	3	2.2	1.7	1.3	1.3
Continuous Flow	1 LPM	2 LPM	3 LPM			
SeQual Eclipse 3 (cont.)	3.7	2.4	1.3			
Respironics SimplyGo (cont.)	2.9	1.6	.9			

Call our customer care agents at **1-877-303-9315** to help you select the best concentrator for your needs.

Note 1: Settings 1-6 are not standardized across manufacturers  
 Note 2: Settings 1-6 on Pulse Dose are not equivalent to 1-6 LPM on Continuous flow