

# Home Oxygen Concentrators

## 2 to 10 Liter Home Concentrators

Unit Name	Liter Flow	Oxygen Concentration	Weight	Dimensions	Noise Level	Power Consumption	Warranty
<b>Respironics SimplyFlo</b> 	Up to 2 LPM	87% to 96%	9 pounds	11.5" Wide 10" High 6" Deep	50 decibals at 2 LPM	120W max	3 Years
<b>AirSep VisionAire</b> 	1- 5 LPM Continuous Flow	87% to 95%	30 pounds	14" Wide 20" High 11.5" Deep	40 decibels	290 Watts	3 Years
<b>Respironics EverFlo</b> 	0.5 - 5 LPM	93% (+/- 3%) @ 5 LPM	31 pounds	15" Wide 23" High 9.5" Deep	40 decibels	350 watts	5 Years
<b>Invacare Perfecto</b> 	0.5 - 5 LPM	87% to 95.6%	43 pounds	13" Wide 23" High 11.5" Deep	43 decibels	280 watts	3 Years
<b>AirSep Intensity</b> 	2 to 10 LPM	2-9 LPM at 89% to 95% 10 LPM at 87% to 93%	58 pounds	16.5" Wide 27.5" High 14.5" Deep	45 decibels	590 watts	3 Years
<b>Invacare Platinum 10</b> 	2 to 10 LPM	87% to 95.6%	53 pounds	18.4" Wide 26.4" High 14.4" Deep	58 decibels	585 watts	3 Years
<b>Respironics M10</b> 	1 to 10 LPM	88% to 96%	53 pounds	19" Wide 27" High 13" Deep	53 decibels	n/a	3 Years